

Protocol for amino allyl labeling

Khan Lab (February, 2005)

RT Reaction:

After TRIzol purification and RNeasy, RNA is re-suspend in H₂O. Take 1 µl for RNA concentration measurement, and make sure the total amount of RNA for labeling is > 15 µg, usually 20 µg. For amplified RNA, use 2 µg amplified RNA per labeling per slide.

Concentrate RNA to 23.2µl with Microcon-30.

	Total RNA	Amplified RNA
RNA	23.2µl	22.2 µl
Anchored primer (1µg/µl)	1 µl	-
(or Random Hexamer for amplified RNA, 3µg/µl)	-	2 µl
<hr/>		
Total volume	24.2µl	24.2µl

70°C for 5 minutes and cool to 42°C

RNase Inhibitor RNAsin	1 µl
5X First Strand Buffer	8 µl
0.1 M DTT	4 µl
50X aa-dNTP mix	0.8 µl
SSII Reverse Transcriptase	2 µl
<hr/>	
Total volume	40 µl

42°C for 30 minutes

Additional 1 µl SSII RT enzyme

42°C for 30 minutes and cool to R/T

Labeling Reaction:

1. To each labeling reaction, add 5 µl of 500mM EDTA (pH8.0), mix well before adding NaOH
2. Add 10 µl 1 M NaOH, mix well
3. 65°C for 20 minutes to hydrolyze RNA, then cool to R/T
4. Add 10 µl 1 M HCl to neutralize NaOH
5. Probe purification with Qiagen PCR spin column protocol.
 - a. Add 350 µl PB buffer to RT reaction and transfer to QIAquick column.
 - b. Spin at maximum speed for 1 minute.
 - c. Wash with 750 µl phosphate wash buffer (5 mM KPO₄, pH 8.0 + 80% EtOH – use this home-made buffer to avoid Tris).
 - d. Spin at maximum speed for 1 minute.
 - e. Repeat wash steps c and d.
 - f. Empty collection tube and spin at maximum speed for 1 minute.

- g. Transfer column to a fresh tube, Elute with 30 μ l elution buffer (4 mM KPO₄, pH 8.5) with 1 minute incubation followed by spin at maximum speed for 1 minute.
 - h. Repeat elution step g.
 - i. Speedvac dry (~40 minutes)
6. Cy Dye coupling
- a. Re-suspend the cDNA in 4.5 μ l 0.1 M carbonate buffer, pH 9.0 (freshly made, every two weeks, store at room temperature).
 - b. Add 4.5 μ l NHS-Cy dyes (re-suspended in DMSO. AVOID ANY MOISTURE!!)
 - c. Incubate the reactions in dark for 1 hour in room temperature.

Purification

1. Add 35 μ l 100 mM NaOAc pH5.2 to the reaction.
2. Add 250 μ l PB buffer and follow QIAquick PCR purification protocol.
3. Elute 2X in 30 μ l H₂O.
At this step, 1 μ l of sample can be taken for dye incorporation analysis using NanoDrop Cy3/Cy5 program.
4. Mix the probes of Cy3 and Cy5 in equal amount. Speedvac dry (~40 minutes).

Pre-hybridization:

1. Prepare pre-hybridization buffer (5xSSC, 0.1% SDS and 1% bovine serum albumin, BSA Sigma Cat# A-9418)
2. Place slides to be analyzed into a staining dish, fill with pre-hybridization buffer, and incubate at 42 °C for 45 minutes
3. Wash the slides in a Wheaton slide rack at room temperature with DIH₂O with agitation for 2 min.
4. Dip the slides in room temperature isopropanol and dry by centrifugation at 700 rpm for 5 min. Inspect slides to make sure there is no protein marks on slides. If there is, repeat step 3 and 4.

Slides should be used immediately following pre-hybridization, at least within one hour.

Hybridization:

1. Re-dissolve the combined probes in 27 μ l 1X hybridization solution (240 μ l H₂O, 250 μ l 20XSSC, 10 μ l 10%SDS, 500 μ l Formamide = 1ml) and combine with COT1-DNA and Poly(A)-DNA to block nonspecific hybridization

Cy5/Cy3 labeled probes	27 μ l
Poly dA (8mg/ml)	1.3 μ l
Yeast tRNA (4mg/ml)	1.3 μ l
CoT-1 DNA (10mg/ml, concentrated)	1.3 μ l
Total	30.9 μ l

2. Denature at 95 °C for 2 minutes
3. Apply 14 μ l 5X SSC onto each modified hybridization chamber to prevent drying.

4. Apply probes to the cover slide that is big enough for the printed area (e.g. 60mm cover slip for 22K array). Then put the pre-hybridized microarray slide on the top of the cover slip, invert, and hybridize in a humidified chamber at 42 °C for 16-24 hours.

Wash:

Place slide into a staining dish with Wash #1 at R/T until cover slip falls off.

Wash #1:	at R/T in 0.2% SDS + 1 X SSC until cover slip falls off
Wash #1:	4 minutes at 42 °C in 0.2% SDS + 1 X SSC
Wash #2:	4 minutes at R/T in 0.2% SDS + 0.1 X SSC
Wash #3:	4 minutes at R/T in 0.06 X SSC

Spin immediately (700 rpm for 5 minutes at R/T)

Appendix:

Master mix for labeling reaction (4):

5X first strand buffer	8 μ l X 4.5 = 36 μ l
0.1 M DTT	4 μ l X 4.5 = 18 μ l
50X aa-dNTP mix	0.8 μ l X 4.5 = 3.6 μ l
RNAse inhibitor RNAsin	1 μ l X 4.5 = 4.5 μ l
Superscript II RTase	2 μ l X 4.5 = 9 μ l
<hr/>	
Total:	15.8 μ l/reaction

Calculation for each 1000ml washing buffer:

0.2% SDS:	20 ml 10%SDS in 1000 ml.
1X SSC:	50 ml 20XSSC in 1000 ml.
0.1XSSC:	5 ml 20XSSC in 1000 ml.
0.06XSSC:	3 ml 20XSSC in 1000 ml.

Solutions:

100 mM aa-dUTP (Sigma, cat # A0410)

Dissolve 10 mg in 191 μ l 0.1 M KPO₄ buffer, pH 7.5. Quantitate this stock solution by diluting an aliquot 1:5000 in the same buffer and measuring A₂₈₉. Stock concentration in mM = OD_{289nm} X 704.

Make 10 ml 0.1 M KPO₄ buffer, pH7.5

1. Prepare 2 solutions: 1 M KH₂PO₄ and 1M K₂HPO₄.

2. Mix:

0.802 ml 1M K₂HPO₄

0.198 ml 1 M KH₂PO₄

3. Bring to 10 ml with H₂O. The pH of this KPO₄ buffer should be at 7.5.

50X aa-dNTP mix

Nucleotide/dye ratio with the ratio of aa-dUTP to dTTP. We use a 2:3 ratio. Prep are a 50X stock:

Start Concentration	Volume	Final Concentration
100 mM dATP	25 µl	25 mM dATP
100 mM dCTP	25 µl	25 mM dCTP
100 mM dGTP	25 µl	25 mM dGTP
100 mM dTTP	15 µl	15 mM dTTP
100 mM dUTP	10 µl	10 mM aa-dUTP

Phosphate wash buffer

1. Prepare 2 solutions: 1 M KH₂PO₄ and 1M K₂HPO₄.

2. Make 1 M KPO₄ (the pH should be 8.5-8.7) mix:

9.5 ml 1M K₂HPO₄

0.5 ml 1 M KH₂PO₄

3. Make 100 ml phosphate washing buffer, Mix in the following order:

15.25 ml H₂O

84.25 ml 95% EtOH

0.5 ml 1 M KPO₄, pH 8.5

(This solution will be slightly cloudy)

Phosphate elution buffer

Dilute 1 M KPO₄ pH 8.5 to 4 mM for Phosphate elution buffer:

1 M KPO₄ pH 8.5 0.4 ml

H₂O 99.6ml

Total 100ml

Carbonate buffer 0.1 M Na₂CO₃, pH 9.0:

Na₂CO₃ 0.27g

H₂O 20ml

Adjust pH to 9.0 with 6 N HCl

Total with H₂O 25 ml

NHS-Cy dyes

NHS-Cy3: AmershamPharmacia Cat# PA23001

NHS-Cy5: AmershamPharmacia Cat#PA25001

1. Resuspend one tube in 73 ul DMSO.
2. Use immediately, or aliquot 4.5 ul into 0.5 ml tubes and store at -80°C. Avoid moistures!!!